

## PCR Plastic Consumables

The Cell Projects range of PCR plastics and consumables are manufactured to strict quality and design criteria all starting with moulds and tools designed and produced by us.

- ✓ *Clean room produced to reduce particle and biological contamination.*
- ✓ *Stringent testing to check integrity, performance, specification and leak free use.*
- ✓ *Bio-Testing to ensure products are free of biological contamination.*
- ✓ *Virgin Plastic for clarity, flexibility and ease of use.*



## 96 well PCR Plates

### Non-Skirted 96 well PCR Plates

Embossed text NS-96-CC

Black printed text NS-96-CB

- ✓ *Raised rim for optimal plate sealing compatible with most cap strips, foils and adhesive film.*
- ✓ *Thin walled for optimal heat transfer.*
- ✓ *Flat and rigid platform to enhance the PCR process.*
- ✓ *Cut off corner for orientation.*
- ✓ *Universal design for 0.2ml Thermal Cyclers.*
- ✓ *Batch tested for physical and biological contamination.*
- ✓ *Colour choices available.*

Non-skirted, unprinted  
NS-96-CC



Non-skirted, black printed  
NS-96-CB



### Full Skirted 96 well PCR Plate FS-96-CC

- ✓ *Same benefits as non skirted plates.*
- ✓ *Compatible with robotic and automated systems.*
- ✓ *Low profile design with raised rims for maximal sealing.*
- ✓ *Sealable with caps, film and foil.*
- ✓ *Stackable.*

Full-skirted  
FS-96-CC




## PCR Plate Sealing Accessories

### 8 Well Strip Caps

Flat Cap sealing strips FC-08-CC


Domed Cap sealing strips DC-08-CC

- ✓ Easy push fit with reliable seal.
- ✓ Transparent with labelling positions.
- ✓ Ends labelled and designed for easy orientation.
- ✓ Optimised for compatibility on the PCR plates.
- ✓ Flat caps can be used for QPCR applications.
- ✓ Clean room produced.



Transparent sealing film TSF-PCR-100P/P

- ✓ Suitable for PCR and QPCR.
- ✓ Adhesive can easily be peeled.
- ✓ Ideal for storage and transport.
- ✓ Seal integrity from -40°C to +80°C.
- ✓ Tested for biological contamination.



## Ordering Information

Cat. No.	Description
<b>PCR Plates</b>	
NS-96-CC	Non-skirted 96 well PCR plates packed 5 per bag, 25 to a box
NS-96-CB	Non-skirted 96 well PCR plates with black printed grid packed 5 per bag, 25 to a box
FS-96-CC	Full-skirted 96 well PCR plates packed 5 per bag, 25 to a box
<b>Accessories</b>	
FC-08-CC	Flat cap 8-well sealing strips packed 12 per bag, 300 to a box
DC-08-CC	Domed cap 8-well sealing strips packed 12 per bag, 300 to a box
TSF-PCR-100P/P	Transparent sealing film for PCR plate, 100 sheets