

<b>Cell Projects Part Code Description</b>	<b>FS-96-CC</b> Full Skirted	<b>NS-96-CC</b> Non-skirted	<b>SS-96-CB</b> Semi-Skirted
	96 Well	96 Well	96 Well
<b>Sequencers</b>			
<b>ABi</b>			
3100 Genetic Analyser		<b>X</b>	<b>X</b>
3130 Genetic Analyser		<b>X</b>	<b>X</b>
310 Genetic Analyser		<b>X</b>	<b>X</b>
3500 Genetic Analyser			
3500			<b>X</b>
3500 XL			<b>X</b>
3700 DNA Analyser		<b>X</b>	<b>X</b>
3730/3730XL DNA Analyser		<b>X</b>	<b>X</b>
<b>Amersham</b>			
MegaBACE™ 500	<b>X</b>		
MegaBACE™ 1000 mark 2	<b>X*</b>		
<b>MJ Research</b>			
BaseStation™	<b>X</b>		
<b>Transgenomic</b>			
Wave	<b>X</b>		

\*Compatible with MegaBACE™ bought after July 2000